



# USB-C HANDLE

Rechargeable Battery Handle with rheotronic®

Industry Leading 10 hours Working Time

The Riester Type C, USB-C rechargeable battery handle with rheotronic® is quickly and easily charged by convenient USB-C charging cable in 3.5 hours. This solid metal handle is made from solid black metal and equipped with a practical bayonet attachment mechanism for quickly interchanging instrument heads.



# All Day Power



Designed for the existing family of Riester Products



ri-derma®  
DERMATOSCOPE



ri-scope®L  
OTOSCOPE



ri-scope®L  
OPHTHALMOSCOPE

A versatile USB-C charging handle providing you long operating time, interchangeable with a wide range of Riester instrument heads including ri-scope® and ri-derma® 2.5 V and 3.5 V

- › Higher capacity, with a 10 hour operating time
- › Battery status indicator visible on handle (bicolor LED indicates fully charged, charging, or battery needs to be charged)
- › Wide range charger (100 - 240 V)
- › Universal USB type C charging technology
- › Fully charged in 3.5 hours



## ORDERING INFORMATION

|                   |  |
|-------------------|--|
| Art.No. 10704.003 | USB Battery Handle Typ C, 3.5 V with rechargeable Li-ionen Battery and Charger, US |
| Art.No. 10704.002 | USB Battery Handle Typ C, 3.5 V with rechargeable Li-ionen Battery and Charger, UK |
| Art.No. 10704.001 | USB Battery Handle Typ C, 3.5 V with rechargeable Li-ionen Battery and Charger, AU |
| Art.No. 10704     | USB Battery Handle Typ C, 3.5 V with rechargeable Li-ionen Battery and Charger, EU |

For the complete line of Riester Products, please visit [Riester.de](https://www.riester.de)